



EL12M

ELECTRIC DYNAMIC TILT TABLE



Dynamic Tilt Table is available in a range of configurations to optimise patient care and outcomes, for functional rehabilitation of handicapped or disabled patients and in cardiac assessments.



Technique Data

- Two electric motor for table Hi-Lo & Tilt function
- Height adjustable: 67cm - 103cm
- Angle of tilting: 0° to 80°
- Angle of headrest: 0° to 70°
- Table length and width: 249cm x 75cm
- Max loading weight: 175KGS to 180KGS
- Certificate: ISO13485/CE /FDA/93 42 EEC
- Wheels: 125mm castors for safety and mobility



V2U HEALTHCARE PTE
Address : 22 Sin Ming Lane #07-88 MidView City,
Singapore 573969

Telephone: +65 65700055
Fax: +65 65703394

Email: business@v2uhealth.com

www.v2uhealth.com