

PHYSIOFIT WATERSHIELD

Cast & Bandage Protector

Ideal for protection of Casts, Bandages, Rashes, Abrasions, Burns, Cuts, Wounds, IV & PICC Lines



**ARM
WATERSHIELD**



**LEG
WATERSHIELD**

FEATURES

- Excellent watertight seal
- Easy one hand application
- Non-Latex reusable design
- Adult and paediatric size available

Model	Description	Loop (OD/ID)	Dimension L X W in cm
330130	Small Arm	15.5/3.5 cm	34.5 x 24 cm
330230	Medium Arm	15.5/3.5 cm	48.5 x 24 cm
330330	Large Arm	15.5/3.5 cm	53.5 x 38 cm
330630	Small Leg	15.5/3.5 cm	32 x 28.5 cm
330530	Medium Leg	15.5/5 cm	44 x 34.5 cm
330430	Large Leg	15.5/5 cm	74.5 x 44 cm

INSTRUCTIONS FOR USE

- Remove cast and bandage protector from the box.
- Pull cast abs bandage over injured part until area is covered.
- Position cast and bandage seal above afflicted area.

WARRANTY: Watershield is latex free and guaranteed to last the life of the cast, bandage, wound, IV and PICC lines.

CONTACT

22 Sin Ming Lane #07-88, Mid View City, Singapore 573969
Tel: +65 65700055 | Email: business@v2uhealth.com

www.v2uhealth.com

