



V2U Duo Stim

V2U Duo Stim is our newest development of full digital stimulator with 32 stimulation modes controlled by press buttons for pain relief and muscle training. There are 5 adjustable TENS modes, 3 adjustable EMS modes and 24 professional presets of TENS and EMS available. This offers quick and effective choice of programs for all user



The V2U Duo Stim TENS/EMS is a fully digital battery operated pulse generator that sends electrical impulses to the nerves and underlying muscle groups. This unit is a combination stimulator of TENS and EMS which can be used for pain relief and muscle stimulation. The device is provided with two controllable output channels, each independent of the other. A pair of electrodes can be connected to each output channel. The intensity level and settings are controlled by press buttons



FEATURES

- Max. 100 mA output, 99 steps of intensity control
- Key Lock Facility
- Large LCD Display
- Pause Function
- Timer Built in for Treatment Compliance
- Automatically shut off in 5 minutes if not used
- 17 Modes of Tens, 15 Modes of EMS, 32 total Modes
- Easy to use Press Buttons Control
- Sturdy Metal Belt Clip
- Low Battery Indicator
- Individual Adjustment of Intensity level
- Treatment Time Countdown

V2U Healthcare Pte Ltd

22 Sin Ming Lane #07-88
MidView City, Singapore 573969
(+65) 65 700055
www.v2uhealth.com