





DYSPHAGIA ELECTRODES

The V2U Healthcare line of dysphagia electrodes are high-quality, low-impedance electrodes that are compatible with existing NMES devices for the treatment of dysphagia. These electrodes offer a cost-effective alternative to your healthcare facility.


PRODUCT FEATURES

- Silver (Ag) coated carbon film for reduced impedance, greater conductivity and uniform current dispersion.
- High-quality, dual layer of gel, formulated for maximum adhesion.
- Gel does not separate from carbon.
- Excellent local customer support.
- Prompt turnaround time on orders.
- Competitive pricing.
- Compatible with existing NMES dysphagia units.
- Cost-effective, clinically acceptable replacement dysphagia electrode.
- Available in our Ultra-Flexible, Blue Foam topcoat.
- Made in the U.S.A.
- Longer snap posts for secure connection from lead wires to electrodes.

V2U Healthcare ELECTRODE CONSTRUCTION



Item Number: 20165-1
(2/pk)



Made in the USA



T: +65 65700055
F: +65 65703394
business@v2uhealth.com