



THERABATH

Indications:

- Arthritis
- Fibromyalgia
- Post immobilisation joint stiffness
- Soft Tissues stiffness due to adhesion formation
- Chronic pain due to sprain and strain
- Chronic soft tissue contractures
- Eczema



Why Therabath?

Paraffin Wax Bath Therapy is a **drug free and non-invasive** method to **relieve pain**.

Scientific Facts:

Paraffin wax absorbs and retains a great amount of heat. Once applied on the body parts, heat is transferred to the core of the affected area as paraffin goes through a phase change.

As it melts, the paraffin becomes a liquid and can retain more heat. Then, when a hand, foot, or other area is dipped into the heated bath, the phase of the paraffin that surrounds the dipped area quickly changes into a solid. The heat that is given off is called the heat of fusion.

Easy to apply on hand

Liquid Wax



After Drying



Pliable ball of wax to effectively exercise the hand during rehabilitation



Treatment of dry and scaly hands to soften the joints and soft tissues and improve mobility



V2U HEALTHCARE pte ltd.

Address : 22 Sin Ming Lane #07-88 MidView City, Singapore 573969

Email : business@v2uhealth.com | Tel : (+65) 65 700055

Website: www.v2uhealth.com



V2U THERABATH

Model: TB7

Made in: USA

First Model: 1962

Wax Capacity of tank: 9lb

STANDARD ACCESSORIES

**TB7-Paraffin wax bath
with 6 Lbs.
(2.72 Kg wax) refill**

Unit dimensions: 14" X 8" X 8", (36 X 20 X 20 cm)
Unit weight: (full) 15lbs (6.80 kg)
Unit capacity: (full) 9lbs (4.08 kg)
Operating Temp.: 124° - 134° F. (52.5° - 56.7° C)

6 Lbs. (2.72 Kg wax) refill

Unit weight: (full) 15lbs (6.80 kg)

Tank bottom grill

Unit capacity: (full) 9lbs (4.08 kg)

**1 - lid with integrated
handle**

Operating Temp.: 124° - 134° F. (52.5° - 56.7° C)

**1 - power cord
1-Therabath Manual**

Safety feature: It has an 'always safe'
temperature range that ensure the paraffin stays
between 126 and 130 degrees



V2U HEALTHCARE pte ltd.

Address : 22 Sin Ming Lane #07-88 MidView City, Singapore 573969

Email : business@v2uhealth.com | Tel : (+65) 65 700055

Website: www.v2uhealth.com